

Specifications

Dimensions: Full size, from 5'5" tall, folded in bag, L x W x H 34" x 12" x 10" Unfolded 38" x 23,5" height variable

Junior, up to 5'5" tall, in bag 3" shorter. Unfolded 3" shorter and 3" narrower

Shopper (lowest weight version) - slightly reduced dimensions due to smaller wheels. In the Junior version it is perfect for children (6+) and short-statured people

Total weight 35 lbs. w/o battery and seat 24 lbs Shopper 26 lbs (18 lbs)

Choice of 3 li-Ion batteries,

Watt hours (Wh)	Weight	Distance on level and 220 lbs load
158 Wh*	2,6 lbs.	9 mls
274,4 Wh*	4 lbs.	15 mls
420 Wh**	5,8 lbs.	22 mls

* Approved for air travel ** Not approved for air travel

Charger 100 to 240 V Charging time totally discharged battery 5 / 8 / 12 hours

Max. recommended load 320 lbs. plus luggage

Motor: Motor: 24 VDC 200 Watt, incorporated in left rear wheel, forward/reverse

Climbing power with 220 lbs. load 15 % (Shopper 6%)

Speed: 4 mph* Shopper 2,8 mph

*Within the EU the top speed for non-licensed motor vehicles operated on public property is limited to 6 kph. Higher speed setting available on request.

Wheels: front 8" x 2" rear 7,25" x 2,5" all filled with dense foam

Brakes: Each rear wheel with band brake and parking brake

Included in each shipment:

1 TravelScoot with battery of choice and charger, duffel bag, fabric "caddy" pouch (for light items; easy to remove and install), tool kit. Cane holder and basket rail are standard.



TravelScoot

Re-Mobilize Your Life!

*At a sensationally low weight of 35 lbs the **Travelscoot** is worldwide by far the lightest and most compact electric mobility scooter. Ideal for everyone who's feet are not playing along but who is otherwise still in a fairly good shape. The innovative folding mechanism allows, for the first time, a 3-wheeled scooter to fold down to the dimensions of a 2-wheeler.*

Set-up and break-down takes less than a minute. It fits completely into its own carrying bag and can be stored where no mobility scooter has ever fit before!

*The **TravelScoot's** snazzy appearance practically eliminates the perception of being a mobility device for the disabled.*

Strangers, if they notice you at all, are more likely to give you a "thumbs up" or "cool ride".

More details you find on our website

www.travelscoot.com

or visit us:

escooter.de

Adi-Maislinger-Str. 7 D-81373 Munich

Open Mo-Fr 10 a.m. to 6.30 p.m.



www.travelscoot.com

Onlineshop: www.travelscootshop.com

www.escooter.de

Tel.: +49-(0)89-57951906

E-Mail: mail@travelscoot.eu



